

# Treat. Monitor. Adapt.

**With Molecular Precision.**

## **Northstar Select™**

Liquid Biopsy with  
More Treatment Options,  
More Often

## **Northstar Response™**

Precise, Tissue-Free  
Treatment Response  
Monitoring

MOLECULAR  
PRECISION  
AT EVERY STEP

Treat.  
Monitor.  
Adapt.

## ● SELECT

### Liquid Biopsy with More Treatment Options, More Often

A tissue-free, ultra-sensitive, 84-gene genomic alteration panel with MSI status. The panel covers all solid tumors alterations associated with FDA-approved targeted therapies and immunotherapies, and provides clinical trial matching.

- 82 genes for single nucleotide variants (SNVs) / indels
- 19 genes for copy number amplifications
- 5 genes for copy number loss
- 9 genes for fusions

#### **An alternative to tissue**

Tissue samples tend to fail at a higher rate, increase TAT, and may not be representative of the entire tumor / malignancy.

### Precise, Tissue-Free Treatment Response Monitoring


Analyzes >500 genomic loci uniquely methylated in cancer cells. Utilizes methylation signature to estimate tumor fraction. A Tumor Methylation Score™ on each report is used to quantify changes in molecule count, enabling treatment response monitoring. The change in methylated tumor molecules across time points can help understand how the tumor is evolving.

#### **A complement to imaging**

Imaging interpretation can be confounded by immunotherapy pseudo-progression.

## ● RESPONSE






## THE TECHNOLOGY

Single-molecule level quantification powered by patented Quantitative Counting Templates™ (QCT™)

BillionToOne's QCT™ technology combines the sensitivity of ddPCR with the breadth of NGS, enabling precise quantification and ultra-sensitive detection of low variant allele fraction (VAF) actionable targets.



# ● SELECT Treat more actionable mutations with ultra-sensitive genomic profiling.

Druggable mutations are actionable at very low VAF. Uncover more treatment options with an ultra-sensitive assay that offers an extremely low limit of detection (LOD). Studies have shown patients with very low-VAF (<0.2%) mutations responded to therapy<sup>1,2</sup>.

50% of alterations fall below 0.5% VAF and are challenging to detect<sup>3</sup>:




Figure 1. A survey of >1,000 plasma samples found that half of all detected SNVs occurred below ~0.5% VAF ●, and a quarter were below 0.2% VAF ●. Variant detection was generally unreliable and variable between assays for alterations lower than 0.5% VAF.<sup>3</sup>

Northstar Select offers ~2X lower LOD<sup>4</sup>, making it possible to uncover more actionable alterations:




Figure 2. Cited LODs of common assays in the market range from 0.25% to 0.50%<sup>3</sup> ●. Northstar Select's LOD ranges from 0.13% to 0.16%<sup>4</sup> ●. LOD is defined as the lowest concentration of an analyte in a sample that can be consistently detected with ≥95% probability.

1. Abida W, Armenia J, Gopalan A, et al. Prospective Genomic Profiling of Prostate Cancer Across Disease States Reveals Germline and Somatic Alterations That May Affect Clinical Decision Making. JCO Precis Oncol. 2017;2017:PO.17.00029. doi:10.1200/PO.17.00029
2. Jacobs MT, Mohindra NA, Shantzer L, et al. Use of Low-Frequency Driver Mutations Detected by Cell-Free Circulating Tumor DNA to Guide Targeted Therapy in Non-Small-Cell Lung Cancer: A Multicenter Case Series. JCO Precis Oncol. 2018;2:1-10. doi:10.1200/PO.17.00318
3. Deveson IW, Gong B, Lai K, et al. Evaluating the analytical validity of circulating tumor DNA sequencing assays for precision oncology. Nat Biotechnol. 2021;39(9):1115-1128. doi:10.1038/s41587-021-00857-z
4. Internal validation data, June 2023

Alteration Type	Reportable Range	Allele Fraction / Copy Number	Analytical Sensitivity	Analytical Specificity
SNVs / Indels	≥0.01%	>0.17% 0.13 - 0.16% 0.06 - 0.12%	>99% 95.5% 78.4%	>99.9%
CNAs*: Amplifications	≥2.1 copies**	2.125 - 2.160 copies***	95%	
CNAs*: Losses	≤1.9 copies	≥1.77 copies	>95%	
Fusions	≥0.02%	>1%	>99%	
MSIs	n/a	0.07%*** >0.4%	95% 100%	

\* This translates to variant detection at a 0.02 - 0.04% allele fraction even with limited DNA input quantities (i.e. 15 - 30ng of cfDNA).

\*\* 2.5% tumor fraction at 6 copies of amplification. \*\*\* Impacted by biological variability.

Table 1. Northstar Select's performance specifications. Internal data on file, June 2023.



MORE  
TREATMENT  
OPTIONS,  
MORE  
OFTEN

SNVs / Indels 82 genes						CNAs: Amplifications 19 genes			Fusions 9 genes	
AKT1	CCNE1	EZH2	JAK2	NOTCH1	RET	AR	ESR1	PIK3CA	ALK	
AKT2	CD274 (PD-L1)	FANCA	JAK3	NPM1	RHOA	BRAF	FGFR1	RAF1	BRAF	
ALK	CDH1	FBXW7	KIT	NRAS	RIT1	CCNE1	FGFR2	RET	FGFR2	
APC	CDK12	FGFR1	KRAS	NTRK1	ROS1	CD274 (PD-L1)	KIT		FGFR3	
AR	CDK4	FGFR2	MAP2K1 (MEK1)	PALB2	SF3B1	CDK4	KRAS		NTRK1	
ARAF	CDK6	FGFR3	MAP2K2 (MEK2)	PDGFRA	SMAD4	CDK6	MET		NTRK2	
ARID1A	CDKN2A	FGFR4	MET	PIK3CA	SMO	EGFR	MYC		NTRK3	
ATM	CDKN2B	GATA3	MLH1	PMS2	STK11	ERBB2	PDGFRA		RET	
BRAF	CHEK2	GNA11	MPL	PTEN	TERT				ROS1	
BRCA1	CTNNB1	GNAQ	MSH2	PTPN11	TP53					
BRCA2	DDR2	GNAS	MSH6	RAD51C	TSC1					
BRIP1	EGFR	HRAS	MTOR	RAD51D	VHL					
CCND1	ERBB2 (HER2)	IDH1	MYC	RAF1						
CCND2	ESR1	IDH2	NF1	RB1						
						CNAs: Losses 5 genes			Biomarker	
						ATM	CDKN2A		Microsatellite Instability Status (MSI)	
						BRCA1	PTEN			
						BRCA2				

Table 2. Northstar Select's gene list.

# ● RESPONSE

Monitor treatment response precisely and adapt your strategy.

## THE POWER OF METHYLATION

This NGS-based test is designed to detect hundreds of genomic loci uniquely methylated in cancer cells. Northstar Response measures the change in methylated tumor molecules, reflecting tumor fraction changes with high accuracy.

### >10x stronger signal than SNV-based ctDNA monitoring assays.

Distinct from a limited number of single nucleotide variants (SNVs), which are commonly used for ctDNA monitoring, Northstar Response interrogates >500 cancer-associated, hypermethylated genomic loci which represents more than 10x the targets. QCT technology combined with the high number of loci and personalized germline filtering enables highly precise quantification.




Figure 3. Median and interquartile range of Tumor Methylation Scores at each tumor fraction. Sheared gDNA from tumor cells were added to sheared gDNA from healthy cells from the same subject at different tumor fraction.

## CLINICAL RELEVANCE

### Northstar Response can aid in the detection of progression.




Figure 4. Northstar Response successfully stratifies ICI immunotherapy clinical outcomes. (A) Change in Tumor Methylation Score from pre-treatment to post-treatment 1,  $p < 0.0001$ . (B) RECIST at post-treatment 1,  $p = 0.14$ .

Source: Ye et al. Methylated ctDNA dynamics correspond with clinical tumor load in metastatic lung cancer patients on therapy. Poster presented at: AACR 2023; April 18, 2023; Orlando, FL. Abstract 5588.


### CASE EXAMPLE 01

Detected partial response 191 days before imaging and subsequently detected progression

#### Stage IV Pancreatic Ductal Adenocarcinoma

Chemotherapy

- folfirinix
- followed by folfiri maintenance
- followed by folfox




### CASE EXAMPLE 02

Detected ctDNA increase 84 days before imaging detected progression

#### Stage IV Colorectal Adenocarcinoma

Chemotherapy + Immunotherapy

- fluorouracil / leucovorin / bevacizumab



MONITOR +  
ADAPT WITH  
PRECISION  
GUIDANCE

## Ease of access to providers and patients

- ✓ Whole blood sampling, no tissue required
- ✓ Appropriate for any late stage solid tumor type
- ✓ Informs and monitors response to targeted therapy & immunotherapy
- ✓ Clear reporting with ~2 week TAT from lab receipt

BillionToOne is committed to making Northstar assays accessible to patients. We have a robust financial support program available and a dedicated billing support team to help answer questions and resolve any concerns. For any questions, please call 833-537-1819 or email [support@northstaronc.com](mailto:support@northstaronc.com).

### **BILLIONTO ONE**

**T** 833.537.1819  
**F** 833.874.0918  
**E** [support@northstaronc.com](mailto:support@northstaronc.com)

[northstaronc.com](http://northstaronc.com)

BillionToOne is a precision diagnostics company founded in 2016 with the mission to make molecular diagnostics more accurate, efficient, and accessible for all. The CLIA-licensed laboratory is based in Menlo Park, CA. BillionToOne has commercialized a series of unique cell-free DNA-based assays, including the first and only single-gene non-invasive prenatal test for recessive conditions and red blood cell fetal antigens, and therapy selection and response monitoring liquid biopsy assays.

CLIA: 05D2275351.  
Northstar Select™, Northstar Response™, and Tumor Methylation Score™ are trademarks of © 2023 BillionToOne, Inc. All rights reserved.

